



OIS WINSTATION ESSENTIAL™

THE OIS WINSTATION ESSENTIAL™ DIGITAL IMAGING SYSTEM IS DESIGNED TO MEET THE NEEDS OF GENERAL OPHTHALMOLOGISTS, WHO NEED A GOOD QUALITY IMAGE, FOR PRIMARY DOCUMENTATION PURPOSES AT A GREAT VALUE. THIS SYSTEM HAS A SIMPLE INTERFACE THAT CONNECTS EASILY TO YOUR FUNDUS CAMERA, USING A 3.2 MEGAPIXEL SENSOR; THIS SYSTEM PROVIDES THE BEST IMAGE QUALITY IN ITS CLASS OF SYSTEMS ON THE MARKET.

Precise Clinical Information at a Great Value





Highlights

- High-resolution 3.2 megapixel sensor (2048X1536 actual image size)
- Instant FA and color fundus imaging in a single, compact digital camera back
- Seamless integration with Ophthalmology Office™—EMR

SYSTEM FEATURES

- AutoMontage™
- AutoImport™
- Contrast, Sharpen, and Zoom
- Negative
- Flip image for laser treatment
- PDT planning
- Automatic image overlay
- Multiple Image Composite
- Image measurement
- Macro—batch processing
- Multiple visit compare
- Annotation
- WebStation™ upgrade available for on-line storage and remote access
- DICOM/VISTA/HL7-compliant

*Some features are optional upgrades



221 Lathrop Way, Suite I
Sacramento, California 95815

Tel: 800.338.8436

Fax: 916.646.0207

www: oisi.com

Specifications**

Computer	Mid tower Intel® Pentium® processor 1-GB RAM or more Dual 250-GB hard drives or larger 3.5" 1.44-MB floppy drive Enhanced keyboard Two-button scroll mouse UPS battery backup
Monitor	19" Flat panel high-resolution color
Archiving	16X/48X DVD-recorder and CD-recorder
Networking and Telemedicine	100/1000 Mb Ethernet 56K V92 modem Remote access software
Image Capture	3.2 megapixel sensor (2048X1536 actual image size), high-resolution B&W and color camera, optical interface for fundus camera, digital sync controller, and interface cabling. Optional 12-Bit CCD for ICG capabilities.
Software	Microsoft Windows® OIS WinStation™ PC Anywhere® Norton AntiVirus®
Dimensions	21"x43" footprint 100V/50-60 Hz–240V/50-60 Hz
Warranty	1 year parts and labor (extended warranty available)
Accessories	Equipment table Equipment cart Printer cart Printer Stereo viewer

*All brand names are the property of their respective owners.

**OIS strives to continuously update and improve upon our technology. Consequently, product specifications may change without notice. Image resolution will vary depending on fundus camera and video adapter in use.