



WINSTATION SOFTWARE

OIS WINSTATION SOFTWARE HAS BEEN DESIGNED TO BE THE MOST COMPLETE IMAGE CAPTURE AND REVIEW TOOL AVAILABLE. OIS' EXTENSIVE BACKGROUND IN OPHTHALMIC PHOTOGRAPHY TRANSLATES TO AN INTUITIVE, USER FRIENDLY INTERFACE. IMAGE ENHANCEMENT FUNCTIONS ARE EASILY ACCESSED BY PULL-DOWN MENUS OR BY TOOL BAR ICONS. MULTIPLE VISITS CAN EASILY BE COMPARED ON A SINGLE SCREEN.

The Most Advanced Imaging Technology Available




1 OIS WinStation Software seamlessly integrates with most fundus cameras. This integration can enable filter control from your computer with only a click of the mouse. This also enables laterality and other data to be embedded in the image as it is taken.



2 Hyper Zoom Window automatically magnifies a customizable region of the image after every shot. This lets you see the detail and focus of an image without any user input, allowing you to critically and efficiently evaluate as you shoot.




3 Easily view recent frames and review in the capture module, enabling you to edit a record without going to the review module.



AutoMontage™

Retinal images seamlessly merged into one color-corrected view, automatically. AutoMontage gives you a wide field view of the retina and the ability to magnify the most subtle retinal details.



Stereo Viewing

Easily advance through stereo pair images of the optic nerve head for 3D viewing with an external stereo viewer.



Multiple Image Composite

Quickly and easily create multiple image composites enabling you to view multiple procedure types, at the same time.



221 Lathrop Way, Suite I
Sacramento, California 95815
Tel: 800.338.8436
Fax: 916.646.0207
www: oisi.com